ACTI^{fit}®

Meniscal scaffold

INSTRUCTIONS FOR USE

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BRUKSANVISNING

BRUKSANVISNING

BRUGSANVISNING

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English

Intended Use

The Actifit meniscal scaffold (Actifit), meniscal ruler accessory and meniscal ruler guide accessory are intended for use in the treatment of irreparable meniscus tears and meniscal tissue loss.

Criteria for Use

- Irreparable, painful and/or dysfunctional meniscal tissue loss requiring partial meniscectomy.
- Surface defect over the entire circumference of the meniscus to ensure stable fixation.
- Skeletally mature man or women, aged 16 to 50 years.
- Well aligned, stable knee joint or ACL deficiencies corrected within 12 weeks after implantation of Actifit.
- ICRS classification ≤ 3.
- No more than 3 prior meniscal surgeries on involved knee.
- Patients understanding and committed to follow the post-operative rehabilitation programme.

Contraindications

- Concomitant PCL insufficiency of the involved knee.
- ICRS classification > 3.
- Uncorrected malformations, axial malalignment in the lower extremity.
- Documented allergy to any of the product components.
- Systemic or local infection.
- Systemic administration of any type of corticosteroid, antineoplastics, immunostimulating or immunosuppressive agent within 30 days before or after implantation of Actifit.
- Evidence of osteonecrosis in the involved knee.
- Medical history that is positive, for, but not restricted to, the following diseases:
 - Rheumatoid arthritis
 - Relapsing polychondritis
 - Severe degenerative osteoarthritis
 - Inflammatory arthritis
- Overweight obese (defined as patients Body Mass Index >35 kg/m²).
- General neurological abnormalities or neurological conditions that tend to preempt the patient's mental ability or willingness to fulfil the requirements of the recommended rehabilitation programme.

Description of the device and the contents of the packaging

Actifit is a single use component biodegradable meniscal scaffold which has been designed to support the body's own physiological pathways of tissue repair by providing a three-dimensional matrix for vascular ingrowth when implanted and attached to the vascularized portion of the meniscus.

Actifit is available in two shapes; medial and lateral. The scaffold is accompanied by two single use accessories, the meniscal ruler and the meniscal ruler guide, both made of a polyamide blend. Both are provided STERILE (Ethylene Oxide (EO) sterilized).

The expected benefits for the patient are pain relief and improved quality of life. Actifit will allow patients with dysfunctional painful meniscus lesions to experience recovery from the dysfunctional meniscus with restoration of its load bearing and shock absorbing functions.

Precautions

- Verify that the Actifit packaging and products are not damaged or otherwise impaired.
- Ensure that the intended knee is being implanted.
- Verify that the labeling on the Actifit packaging that the correct device (medial or lateral) is being used for the knee to be implanted.
- Do not use the Actifit meniscal ruler without a meniscal ruler guide.
- Open the inner sterile pouch slowly and with care to ensure safe removal of the tray containing the meniscal scaffold.
- Do not clean the meniscal ruler or the meniscal ruler guide with alcohol.
- Do not use sutures other than those recommended in these Instructions for Use.
- Safety for use during pregnancy or delivery has not been established.
- Safety and performance of Actifit have not been established for patients under 16 years and above 50 years of age.
- Do not apply excessive force to the meniscal ruler and/or the meniscal ruler guide. Reposition the patient's lower limb if the meniscal ruler and/or the meniscal ruler guide cannot move freely.

Warnings

- Careful handling is required to avoid damage to the patient's tissue and articular cartilage and the scaffold, the meniscal ruler and/or the meniscal ruler guide.
- Actifit must be stored in the original sterile packaging, unopened.
- Actifit, the meniscal ruler and the meniscal ruler guide are designed for single use only. Actifit has been sterilized using Ethylene Oxide (EO) gas.
- Do not re-sterilize.
- Respect sterility rules during handling.
- Do not implant Actifit if the sterile package has been pierced, damaged or altered or if the expiration date (Use by Date), indicated on the barcode or on the expiry label, has past.
- The products (including any residue, unused or unsterile Actifit device material, the meniscal ruler and meniscal ruler guide) should be disposed of immediately in accordance with prevailing hospital procedures for hazardous material.
- Any opened or damaged packaging should be returned to Orteq Sports Medicine Ltd.
- Sutures when not placed correctly may cause pain after the surgery.
- Damaged, pierced or altered packaging, use after the expiration date, re-use and/or re-sterilization of Actifit, the meniscal ruler or the meniscal ruler guide may cause device defects or degradation resulting in Actifit not performing as intended.
- Non-compliance with the specific instructions for cutting Actifit can lead to suboptimal product performance and/or fixation of the device.
- Sutures should be stable however overly tight suturing may lead to collapse of the porous structure and may subsequently result in suboptimal tissue ingrowth.

Surgical Procedure

Only orthopaedic surgeons experienced in meniscal suturing and who have received adequate training on handling the device should use the Actifit scaffold, the meniscal ruler and the meniscal ruler guide. The device is intended for partial meniscus defects extending to the vascularized zone with sufficient anterior and posterior horn tissue as well as meniscus rim tissue to ensure secure fixation. Standard arthroscopic surgery procedure for meniscectomies and standard equipment are applicable for the implantation of Actifit.

Surgical debridement is performed to carefully remove any damaged or pathological tissue and to straighten edges of the defect in the damaged portion of the meniscus to ensure an optimal fit and good connection with the vascularized part of the native meniscus. Care should be taken at all times not to remove more tissue than called for. The accompanying meniscal ruler in combination with the meniscal ruler guide is used to assess the size of the defect post debridement and subsequently estimate the appropriate sizing of Actifit.

The meniscal ruler is inserted through the meniscal ruler guide from the proximal to the distal end of the meniscal ruler guide. Depending on the clinical situation the ruler markings can either be on the inside or the outside of the curve of the meniscal ruler guide.

The meniscal ruler guide with the meniscal ruler inserted are inserted into the appropriate arthroscopic portal in the knee. Sufficient portal width should be ensured to enable easy access. Position the tip of the meniscal ruler against the posterior side of the meniscus defect, the ruler markings should be clearly seen on the arthroscopic tower screen. When in position, the meniscal ruler is deployed by pushing the meniscal ruler through the meniscal ruler guide. Determine the meniscus defect length. Once the length of the meniscus defect is assessed, the meniscal ruler and the meniscal ruler guide are removed. Use the same meniscal ruler to determine the appropriate length of the scaffold and therefore where it should be cut.

Actifit should be cut with a scalpel and in such a place and manner as to ensure that sterility is maintained at all times. Align the meniscal ruler tip at the proximal side of the scaffold and curve the meniscal ruler along the outside circumference of the scaffold. For the purpose of achieving a good fit into the defect the length of Actifit should be oversized by approximately 3 mm for small defects (<3 cm) and approximately 5 mm for large defects (≥3 cm). It is recommended that the anterior edge should be cut in an angle of approximately 45 degrees to facilitate suturing and sutured last.

The tailored Actifit device is then brought into the knee joint through the anteromedial or anterolateral portal using a non-cannulated tissue tension grasper such as the Acuflex Grasper Tissue Tensioner™ (Smith & Nephew) or any other atraumatic

grasper. Sufficient width of the insertion portal should be ensured.

Finally, Actifit should be fixed to the native meniscus rim and horns using standard meniscus suture techniques, such as outside-in or inside-out until secure fixation has been achieved. Standard commercially available size 2.0 non-resorbable sutures are recommended, such as polyester or polypropylene and braided or monofil sutures. In addition, all-inside meniscal repair devices such as FAST FIX™ / Ultra FAST-FIX™, FAST FIX 360 or any other suitable all-inside device may be used.

Home or Occupational Environment

Patients must follow the rehabilitation guidelines. Failure to do so may result in knee injury.

Adverse events

Possible adverse events include but are not limited to the following:

- Transient or permanent nerve damage due to pressure or haematoma.
- Complications associated with arthroscopic suturing technique.
- Bleeding and/or wound irritation and/or slow wound healing.
- Pain.
- Swelling.
- Restricted freedom of movement.
- Any complication associated with arthroscopic surgery.

Principal risks related to the use of the device

- Dislocation of the device after implantation.
- Absence of function; the chance of a patient requiring removal of the Actifit® following implantation of the device is approximately 2%.
- Allergic reactions to the material and/or its degradation products.

Disclaimer of warranty

Although Actifit and the meniscal ruler and meniscal ruler guide hereafter collectively referred to as "Product", have been carefully designed, manufactured and tested prior to use, the Product may fail to perform its intended function satisfactorily for a variety of reasons. Warnings contained in the Product labeling and these Instructions for Use provide more detailed information and are considered an integral part of this disclaimer. Orteq Sports Medicine Ltd therefore, disclaim all warranties, both express and implied, with respect to the Product. Orteq Sports Medicine Ltd shall not be liable for any incidental or consequential damages caused by any use, defect or failure of the Product, whether the claim is based on warranty, contract, tort or otherwise. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the disclaimer of warranty shall not be affected and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

Actifit is made of an aliphatic polyurethane consisting of poly (ε-caprolactone), butanediisocyanate (BDI) and butanediol (BDO). The end degradation products of the Actifit polyurethane are hydroxy hexanoic acid (HHA), BDO and butanediamine (BDA, also known as putrescine).

The meniscal ruler and meniscal ruler guide are made of polyamide and are non-degradable and non-recyclable.

The Product and packaging are latex-free.

© 2018 Orteq Sports Medicine Ltd., All rights reserved

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© 2018 Orteq Sports Medicine Ltd. All rights reserved